

Product Codes of Most Frequently Distributed Products

Code	Red Blood Cells
E0181	RED BLOOD CELLS CPD/500mL/refg ResLeu:<5E6
E0332	RED BLOOD CELLS CPD>AS1/500mL/refg Irradiated ResLeu:<5E6
E0336	RED BLOOD CELLS CPD>AS1/500mL/refg ResLeu:<5E6
E4526	Apheresis RED BLOOD CELLS ACD-A>AS1/XX/refg Irradiated ResLeu:<5E6
E4527	Apheresis RED BLOOD CELLS ACD-A>AS1/XX/refg Irradiated ResLeu:<5E6 1st container
E4528	Apheresis RED BLOOD CELLS ACD-A>AS1/XX/refg Irradiated ResLeu:<5E6 2nd container
E4531	Apheresis RED BLOOD CELLS ACD-A>AS1/XX/refg ResLeu:<5E6
E4532	Apheresis RED BLOOD CELLS ACD-A>AS1/XX/refg ResLeu:<5E6 1st container
E4533	Apheresis RED BLOOD CELLS ACD-A>AS1/XX/refg ResLeu:<5E6 2nd container
E4538	Apheresis RED BLOOD CELLS ACD-A>AS3/XX/refg Irradiated ResLeu:<5E6
E4539	Apheresis RED BLOOD CELLS ACD-A>AS3/XX/refg Irradiated ResLeu:<5E6 1st container
E4540	Apheresis RED BLOOD CELLS ACD-A>AS3/XX/refg Irradiated ResLeu:<5E6 2nd container
E4543	Apheresis RED BLOOD CELLS ACD-A>AS3/XX/refg ResLeu:<5E6
E4544	Apheresis RED BLOOD CELLS ACD-A>AS3/XX/refg ResLeu:<5E6 1st container
E4545	Apheresis RED BLOOD CELLS ACD-A>AS3/XX/refg ResLeu:<5E6 2nd container
E5157	RED BLOOD CELLS CPD>AS1/XX/refg ResLeu:<5E6 LowVol:anticoag not adj
E5122	RED BLOOD CELLS CPD>AS1/XX/refg Irradiated ResLeu:<5E6 LowVol:anticoag not adj
	Auto and Dir RBC Products may include any RBC codes with collection type designation in product suffix.

Code	Plasma Products
E2457	Liquid PLASMA CPD/XX/refg
E2462	Liquid PLASMA CPD/XX/refg Irradiated
E2553	PLASMA CPD/XX/<=-18C Cryo reduced
E2555	PLASMA CPD/XX/<=-18C Frozen <=24h
E7644	Apheresis PLASMA ACD-A/XX/<=-18C RT<=24h frozen<=24h
E7646	Apheresis PLASMA ACD-A/XX/<=-18C RT<=24h frozen<=24h 1st container
E7648	Apheresis PLASMA ACD-A/XX/<=-18C RT<=24h frozen<=24h 2nd container
E7650	Apheresis PLASMA ACD-A/XX/<=-18C RT<=24h frozen<=24h 3rd container
E7607	Apheresis PLASMAJACD-A/XX/<=-18CJRT<=24h frozen<=24hJ4th container

Code	Cryoprecipitate
E5165	CRYOPRECIPITATE None/XX/<=-18C
E3587	POOLED CRYOPRECIPITATE None/XX/<=-18C



Product Codes of Most Frequently Distributed Products

Code	Whole Blood Derived Platelet (Pool)
EA045	POOLED PLATELETS CPD/XX/20-24C ResLeu:<5E6 Bacterial monitoring >=36h
EA553	POOLED PLATELETS CPD/XX/20-24C Irradiated ResLeu:<5E6 Bacterial monitoring >=36h
EA789*	Pooled Platelets CPD/XX/20-24C ResLeu:<5E6 4 units Bacterial monitoring
EA790*	Pooled Platelets CPD/XX/20-24C ResLeu:<5E6 5 units Bacterial monitoring
EA791*	Pooled Platelets CPD/XX/20-24C ResLeu:<5E6 6 units Bacterial monitoring
EA792*	Pooled Platelets CPD/XX/20-24CIIrradiated ResLeu:<5E6 4 units Bacterial monitoring
EA793*	Pooled Platelets CPD/XX/20-24CIIrradiated ResLeu:<5E6 5 units Bacterial monitoring
EA794*	Pooled Platelets CPD/XX/20-24CIIrradiated ResLeu:<5E6 6 units Bacterial monitoring
EA795*	Pooled Platelets CPD/XX/20-24C ResLeu:<5E6 <3E11 plts Bacterial monitoring
EA796*	Pooled Platelets CPD/XX/20-24C Irradiated ResLeu:<5E6 <3E11 plts Bacterial monitoring

*lowa and Missouri facilities ONLY

Code	Single Donor Platelets
E5030	Apheresis PLATELETS ACD-A/XX/20-24C ResLeu:<5E6 Bacterial monitoring
E5031	Apheresis PLATELETS ACD-A/XX/20-24C ResLeu:<5E6 1st container Bacterial monitoring
E5032	Apheresis PLATELETS ACD-A/XX/20-24C ResLeu:<5E6 2nd container Bacterial monitoring
E5033	Apheresis PLATELETS ACD-A/XX/20-24C ResLeu:<5E6 3rd container Bacterial monitoring
E5034	Apheresis PLATELETS ACD-A/XX/20-24C Irradiated ResLeu:<5E6 Bacterial monitoring
E5035	Apheresis PLATELETS ACD-A/XX/20-24C Irradiated ResLeu:<5E6 1st container Bacterial monitoring
E5036	Apheresis PLATELETS ACD-A/XX/20-24C Irradiated ResLeu:<5E6 2nd container Bacterial monitoring
E5037	Apheresis PLATELETS ACD-A/XX/20-24C Irradiated ResLeu:<5E6 3rd container Bacterial monitoring
E5073	Apheresis PLATELETS ACD-A/XX/20-24C Irradiated ResLeu:<5E6 <3E11 plts Bacterial monitorin
E5075	Apheresis PLATELETS ACD-A/XX/20-24C ResLeu:<5E6 <3E11 plts Bacterial monitoring

Code	Pathogen Reduced Platelets
E8331	Apheresis PLATELETS ACD-A/XX/20-24C ResLeu:<5E6 Psoralen-treated
E8332	Apheresis PLATELETS ACD-A/XX/20-24C ResLeu:<5E6 1st container Psoralen-treated
E8333	Apheresis PLATELETS ACD-A/XX/20-24C ResLeu:<5E6 2nd container Psoralen-treated



Product Codes of Most Frequently Distributed Products

Code	Washed RBC Products
E4144	Washed Apheresis RED BLOOD CELLS None/XX/refg Open ResLeu:<5E6
E4145	Washed Apheresis RED BLOOD CELLS None/XX/refg Open Irradiated ResLeu:<5E6
E4560	Washed Apheresis RED BLOOD CELLS None/XX/refg Open Irradiated ResLeu:<5E6 1st container
E4561	Washed Apheresis RED BLOOD CELLS None/XX/refg Open Irradiated ResLeu:<5E6 2nd container
E4564	Washed Apheresis RED BLOOD CELLS None/XX/refg Open ResLeu:<5E6 1st container
E4565	Washed Apheresis RED BLOOD CELLS None/XX/refg Open ResLeu:<5E6 2nd container
E5169	Washed RED BLOOD CELLS None/XX/refg Open ResLeu:<5E6
E5170	Washed RED BLOOD CELLS None/XX/refg Open Irradiated ResLeu:<5E6

Code	Washed Platelet Products
E3541	Washed Apheresis PLATELETS None/XX/20-24C Open Irradiated ResLeu:<5E6
E3543	Washed Apheresis PLATELETS None/XX/20-24C Open Irradiated ResLeu:<5E6 1st container
E3544	Washed Apheresis PLATELETS None/XX/20-24C Open Irradiated ResLeu:<5E6 2nd container
E3545	Washed Apheresis PLATELETS None/XX/20-24C Open Irradiated ResLeu:<5E6 3rd container
E3556	Washed Apheresis PLATELETS None/XX/20-24C Open ResLeu:<5E6
E3558	Washed Apheresis PLATELETS None/XX/20-24C Open ResLeu:<5E6 1st container
E3559	Washed Apheresis PLATELETS None/XX/20-24C Open ResLeu:<5E6 2nd container
E3560	Washed Apheresis PLATELETS None/XX/20-24C Open ResLeu:<5E6 3rd container

Code	Deglycerolyzed RBC Products
E4519	Deglycerolized RED BLOOD CELLS None/XX/refg Open ResLeu:<5E6
E4521	Deglycerolized RED BLOOD CELLS None/XX/refg Open Irradiated ResLeu:<5E6
E4583	Deglycerolized Apheresis RED BLOOD CELLS None/XX/refg Open Irradiated ResLeu:<5E6 1st container
E4584	Deglycerolized Apheresis RED BLOOD CELLS None/XX/refg Open Irradiated ResLeu:<5E6 2nd container
E4588	Deglycerolized Apheresis RED BLOOD CELLS None/XX/refg Open ResLeu:<5E6 1st container
E4589	Deglycerolized Apheresis RED BLOOD CELLS None/XX/refg Open ResLeu:<5E6 2nd container